EN Safety Precautions

Safety Precautions

Please follow the rules outlined in this guide closely to ensure your safety and flawless operation:

- This product is only intended for use on the face.
- **Discard** using the product immediately if any abnormal situation, allergic reaction or discomfort occurs during use.
- This product is a skin care device and is not made for medical purposes. If you are under medical treatment, please consult with your doctor before use.
- This product is not suitable to use on body parts that are embedded with silicone, metal or other foreign objects. Do not use directly on the eyes.
- It is recommended to wait at least 6 months after any facial treatments such as a face lift or botox injections before using this product. If you are unsure, consult with your doctor.
- Thoroughly clean the product after use.
- Do not use any kind of liquid to clean the product. Keep in mind that the surface of the product should only be wiped with a clean dry cloth.
- This product has a special lacquer - when cleaning, do not use any abrasive solution which may cause damage to the lacquer or appearance.
- Replace the exfoliating head every three months.
- Repair, removal or modification of this product is strictly prohibited. Improper handling may result in short circuits, product damage or even fire. Maintenance and warranty services will be void if unauthorized changes have been made to the product. Please contact an authorized dealer to obtain after-sales service if the product has a problem.
- Be cautious of any edges, uneven surfaces, metal parts, accessories or packaging to prevent possible injury or damage. Take note that this product should not be placed in your mouth.
- Do not use this product in locations with high temperatures or high humidity levels (e.g., bathrooms); ideally, the product should be used in a location free of excessive dust and not exposed to direct sunlight.
- Do not insert foreign objects or liquids into the product; this could damage the product.
- Safely from dropping the product or bumping it against other objects. Avoid damage, collision and excessive pressure by always handling the product with care.
- Do not place this product near strong magnetic or electric fields - this may cause the product failure.
- Do not use this product in hazardous locations.
- The product may explode if disposed in a fire.
- If the unit will not be using for an extended period of time, please charge the unit once every three months.
- Please do not place the unit in a damp place like bathroom or sink basin for a long period of time.
- Keep away from water when charging, do not charge near a bathtub, washbasin, or other utensils.

Children

Keep out of reach of children. Do not allow children to play with your product or its accessories. They could hurt themselves or others and could accidentally damage the product. Your product and its accessories may contain small parts that could become detached and create a choking hazard.

Disposal of Electrical & Electronic Equipment

This symbol indicates that all electrical and electronic equipment included with this product should not be treated as household waste. Instead, it should be left at the appropriate collection point for the recycling of electronic equipment. By ensuring this product is disposed of correctly, you are helping to prevent negative consequences for the environment and human health. To recycle this product, please contact your local city office, environmental protection agency, household waste disposal service or where you purchased the product to learn more about appropriate recycling methods.

Battery

It is recommended that you completely use the initial charge up, then recharge the product at least two more times to ensure the best possible performance.

Limited Warranty

If within 12 months of purchase this product or any part thereof is proven to be defective by reason of faulty workmanship or defect in materials, we will repair or replace the same product free of charge of labor or materials on condition that:

1. The original invoice or sales receipt with the date of purchase is presented with the defective unit.
2. Lifetrons does not take any responsibility for lost information of any kind and will not reimburse you for any such loss.
3. The limited warranty is only valid if the original proof of purchase for this product, issued by a Lifetrons authorized dealer, specifying the date of purchase, is presented with the product to be repaired or replaced. Lifetrons reserves the right to refuse warranty service if this information has been removed or changed after the original purchase of the product from the dealer.
4. A rechargeable battery can be charged and discharged over one hundred times. However, it will eventually wear out - this is not a defect and corresponds to normal everyday use. When the operating or standby time is noticeably shorter, it may be a sign that the battery has reached the end of its life span. Please do not make repairs, modifications, or installations without contacting an authorized Lifetrons representative.
5. The warranty does not cover product failures which have been caused by the use of accessories or other peripheral devices which are not Lifetrons branded original accessories intended for use with the product. The warranty does not cover any damage caused by misuse, abuse, neglect or accidental breakage.
6. The warranty does not cover product failures caused by modifications, repairs or disassembly by any person not authorized by the manufacturer.
7. Any defective unit or part shall become property of the manufacturer.

Additional Information

Register your Lifetrons Beauté product online at www.lifetrons.com.

Please visit our website www.lifetrons.com or contact our customer service center at customer.cs@lifetrons.com to learn more about our new products and special offers.

Product specifications are subject to change at the sole discretion and decision of Lifetrons Switzerland AG, Lifetrons International Ltd. or Li-Tek without prior notice. Any devices or accessories not listed under "What's Inside" are in reference to compatibility only and are not included. Product images are for reference purposes only, color may vary from the actual product.

"Lifetrons Beauté", "Lifetrons Beauté Switzerland", "Lifetrons Beauté Professional Skincare", "Lifetrons Switzerland", "Beauté Switzerland", "MyB" are trademarks or registered trademarks owned by Lifetrons Switzerland AG, Lifetrons International Ltd. or Li-Tek. All other brands

and product names listed herein are trademarks or registered trademarks owned by their respective holders and are hereby acknowledged and recognized as such. ©Copyright 2016 by Lifetrons Switzerland AG, Lifetrons International Ltd. and Li-Tek. All rights reserved.

This user guide was originally composed in English and was subsequently translated into other languages. The fidelity of subsequent translations cannot be guaranteed. In case of conflict between the English version and another language version, the English version takes precedence.

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"Lifetrons Beauté", "Lifetrons Beauté Switzerland", "Lifetrons Beauté Professional Skincare", "Lifetrons Switzerland", "Beauté Switzerland", "MyB" are trademarks or registered trademarks owned by Lifetrons Switzerland AG, Lifetrons International Limited and Li-Tek. All other brands and product names listed herein are trademarks or registered trademarks owned by their respective holders and are hereby acknowledged and recognized as such. ©Copyright 2016 by Lifetrons Switzerland AG, Lifetrons International Ltd. and Li-Tek. All rights reserved.

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"Lifetrons Beauté", "Lifetrons Beauté Switzerland", "Lifetrons Beauté Professional Skincare", "Lifetrons Switzerland", "Beauté Switzerland", "MyB" are trademarks or registered trademarks owned by Lifetrons Switzerland AG, Lifetrons International Limited and Li-Tek. All other brands and product names listed herein are trademarks or registered trademarks owned by their respective holders and are hereby acknowledged and recognized as such. ©Copyright 2016 by Lifetrons Switzerland AG, Lifetrons International Limited and Li-Tek. All rights reserved.

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CM-400

User Guide

Thank you for purchasing this Lifetrans® product. This user guide outlines the operating procedures and safety precautions. Please read it carefully to make sure you are using your product properly and retain it for future reference.

How to Use

A. Facial cleaning mode

1. Please remove the make up before using the facial cleaning mode.
2. Press button ① to get start, default facial cleaning function with head ②, at the same time indicator ③ lights up in red, indicator ④ default as gear 1.
3. Press button ⑤ to increase the intensity, press button ⑥ to decrease the intensity (see indicator ③).

Strength Level

1. LED lights up as gear 1
2. LEDs light up as gear 2
3. LEDs light up as gear 3
4. LEDs light up as gear 4
5. LEDs light up as gear 5

4. Select the appropriate strength level for your skin, then start the cleaning process, please see below pictures and follow.
- 1). Apply a foaming cleanser.
- 4.2). Clean the cheeks.
- 4.3). Clean the forehead.
- 4.4). Clean the skin of nose.
- 4.5). Clean the lower part skin of the eyes.

B. EMS with vibration function mode

1. Please remove the make up and clean the skin before using.
2. To ensure proper operation and achieve the maximum effectiveness, use serum with the following procedures.
3. Press button ① to get start, when the unit turns on, press button ⑤ again to change to EMS & vibration mode, a beep sounds indicated, at the same time indicator ③ lights up in blue, indicator ④ default as gear 1 hold the device with touching the conductive pad ⑦, then the micro-current function will be activated, and you will feel a bit electric-pinched, but it's a normal phenomenon.
4. Press button ⑥ to increase the intensity, press button ⑤ to decrease the intensity (see indicator ③).

Strength Level

1. LED lights up as gear 1
2. LEDs light up as gear 2
3. LEDs light up as gear 3
4. LEDs light up as gear 4
5. LEDs light up as gear 5

5. Follow the pictures below, use head ② and follow the texture to do the treatment.
- 5.1). Massage along the face from the lower jaw upward to lobe direction.
- 5.2). Massage along the wing of nose upward to the temples direction.
- 5.3). Massage along the arrow direction.

C. Ultrasonic & ionic negative import function mode

1. Please remove the make up and clean the skin before using.
2. To ensure proper operation and achieve the maximum effectiveness, use serum or facial mask with the following procedures.
3. Press button ① to get start, when the unit turns on, press button ⑤ twice to change to Ultrasonic & ionic negative import function mode, a beep sounds indicated, at the same time indicator ③ lights up in green, indicator ④ default as gear 1 hold the device with touching the conductive pad ⑦, then the negative-ion function will be activated.
4. Press button ⑥ to increase the intensity, press button ⑤ to decrease the intensity (see indicator ③).

Strength Level

1. LED lights up as gear 1
2. LEDs light up as gear 2
3. LEDs light up as gear 3
4. LEDs light up as gear 4
5. LEDs light up as gear 5

5. Follow the pictures below, use head ② and follow the texture to do the treatment.
- 5.1). Massage along the face from the lower jaw upward to lobe direction.
- 5.2). Massage along the wing of nose upward to the temples direction.
- 5.3). Massage along the arrow direction.

D. How to turn the unit off

The unit will automatically turn off after 15 minutes continuous use. A buzzer sounds once per minute in order to trace for treatment time, (total 14 times). If you want to shorten the use of time, press and hold the power button for 3 seconds, there will be a buzzer sounds in 3 times before the unit turns off.

E. Properly charge

1. Piece the unit onto the charge stand, make sure a well connection of part ⑩ and ⑪.
2. Plug the head of ⑫ (Micro USB interface) to head of ⑬ (Micro USB port).
3. Plug the head of ⑭ (USB interface) to a USB port with power supplier (e.g. an AC adaptor with USB port)

F. LED Indicator

- Lower battery:**
- 1). If the vibration intensity is not enough, means the battery is low, please charge the unit.
 - 2). Indicator ③ (in red) with flash in every 5 seconds.
- Battery charging:** Indicator ③ flash per second.
- Battery fully charged:** Indicator ③ is permanently on for total 5.

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