Ultra Photon Light
Red & Blue Light Skin Treatment

EN User Guide

Thank you for purchasing this Lifetrons® product. This user guide outlines the operating procedures and safety precautions. Please read it carefully to make sure you are using your product properly and retain it for future reference. This product is applicable for face only.
② Refer to instruction manual booklet

How to Use

A. Select photon light head

Choose your own color head ① (red head) or ② (blue head).

B. Red photon light mode

1. Please use makeup remover and cleanse the skin.
2. Hold the main unit on the bottom with right hand and pinch the ① by the left hand. Load the unit on the top with the head and hear the pop sound when it fits.
3. Wear the goggles.
4. Press ② (blue LED is always on) DI to start, ① is on, defaults to the last selected working mode.
5. Move the ① from about 5 cm to the face, and move slowly from bottom to top.

C. Blue photon light mode

1. Please use makeup remover and cleanse the skin.
2. The right hand holding the host, left hand pinch ② into the host, pull it out, hear pop sound once been sucked in place.
3. Wear the goggles.
4. Press ③ (blue LED is on) with DI sound to start, ② always on, the default mode is selected.
5. Hold ② away from the facial skin about 5 cm, slowly move it from bottom to top.

D. Treatment time setting

1. In the standby mode, press ②, the first level blue LED will be on. After 1 minute, the device will beep three times consecutively.
2. In the first level position, press ②, the second level blue LED will be on. After 3 minutes, the device will beep for three times consecutively.
3. In the 2nd level position, press ②, the third level blue LED will be on. After 5 minutes, the device will beep for three times consecutively.
4. In the third level position, press ②, the fourth level blue LED will be on. After 12 minutes, the device will beep for three times consecutively.
5. In the fourth position, press ②, the L level blue LED will be on. After 120 minutes, the device will beep for three times consecutively.
6. In the L position, press ② to enter standby mode. After 1 minute, the device will beep for three times consecutively.

E. Pulse light mode

1. In any operating mode, press ③ once (blue LED is on) to enter pulsed mode, ① blinks.
2. Press ③ (blue LED goes off) again for a "DI" sound, you can return to the working mode, ① is on.

F. Shutdown

1. In any operating mode, press ①, it will beep three times consecutively and then turn off.
2. In any operating mode, after setting the timing, the device will beep three times consecutively after closing it.

G. Charging

1. Connect ④ (mini USB plug) and ⑤ (mini USB connector).
2. Connect ⑥ (USB plug) to a powered USB port (for example, an AC adapter with a USB power supply).

H. Battery Indicator

Low Battery: ⑦ (Green LED) blinks
During charging: ⑧ (red LED) is on
Battery full: ⑨ (red LED) goes out

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E. Mc
1. In c
suc
puli
2. Pre
tor
acc

Safety Precautions

Please follow the rules outlined in this guide closely to ensure your safety and flawless operation:

- Be cautious of any edges, uneven surfaces, metal parts, accessories or packaging to prevent possible injury or damage. Take note that this product should not be placed in your mouth.
- Do not modify, repair or dismantle the product. Doing so may result in fire, electric shock, complete breakdown of the product, etc. and will void the warranty.
- Do not use any kind of liquid to clean the product. Keep in mind that the surface of the product should only be wiped with a clean dry cloth.
- Do not use this product in locations with high temperatures or high humidity levels (e.g., bathrooms); ideally, the product should be used in a location free of excessive dust and not exposed to direct sunlight.
- Do not insert foreign objects or liquids into the product; this could damage the product.
- Refrain from dropping the product or bumping it against other objects.
- Do not use this product in hazardous locations.
- Should any abnormal situations or problems occur while using this product, discontinue its use immediately.
- The product may explode if disposed of in a fire.
- Immediately discontinue using this device if an allergic reaction occurs or you feel discomfort.
- This product is not suitable to use on body parts that are embedded with silicone, metal or other foreign objects. Do not use directly on the eyes.
- It is recommended to wait at least 6 months after any facial treatments such as a face lift or Botox injections before using this product. If you are unsure, consult with your doctor.
- This product is a skin care device and is not made for medical purposes. If you are under medical treatment, please consult with your doctor before use.
- Do not place this product near strong magnetic or electric fields - this may cause product failure to occur.
- This product has a special lacquer - when cleaning, do not use any corrosive solution which may cause damage to the lacquer or appearance.

Children

Keep out of reach of children. Do not allow children to play with your product or its accessories. They could hurt themselves or others and could accidentally damage the product. Your product and its accessories may contain small parts that could become detached and create a choking hazard.

Disposal of Electrical and Electronic Equipment

This symbol indicates that all electrical and electronic equipment included with this product should not be treated as household waste. Instead, it should be left at the appropriate collection point for the recycling of electronic equipment. By ensuring this product is disposed of correctly, you are helping to prevent negative consequences for the environment and human health. To recycle this product, please contact your local city office, household waste disposal service or where you purchased the product.

Limited Warranty

If within 12 months of purchase this product or any part thereof is proven to be defective by reason of faulty workmanship or defect in materials, we will repair or replace the same product free of charge of labor or materials on condition that:

1. The original invoice or sales receipt with the date of purchase is presented with the defective unit.
 2. Lifetrons does not take any responsibility for lost information of any kind and will not reimburse you for any such loss.
 3. The limited warranty is only valid if the original proof of purchase for this product, issued by a Lifetrons authorized dealer, specifying the date of purchase, is presented with the product to be repaired or replaced. Lifetrons reserves the right to refuse warranty service if this information has been removed or changed after the original purchase of the product from the dealer.
 4. The warranty does not cover product failures which have been caused by use of accessories or other peripheral devices which are not Lifetrons branded original accessories intended for use with the product. The warranty does not cover any damage caused by misuse, abuse, neglect or accidental breakage.
 5. The warranty does not cover product failures caused by modifications, repairs or disassembly by any person not authorized by the manufacturer.
3. Any defective unit or part shall become property of the manufacturer.

Additional Information

Register your Lifetrons Beauté product online at www.lifetrons.com.

Please visit our website www.lifetrons.com or contact our customer service center at cs@lifetrons.com to learn more about our new products and special offers.

Product specifications are subject to change at the sole discretion and discretion of Lifetrons Switzerland AG, Lifetrons International Ltd. or Li-Tek without prior notice. Any devices or accessories not listed under "What's Inside" are in reference to compatibility only and are not included. Product images are for reference purposes only, color may vary from the actual product.

"Lifetrons" is a trademark or registered trademark owned by Lifetrons Switzerland AG or Lifetrons International Limited. All other brands and product names listed herein are trademarks or registered trademarks owned by their respective holders and are hereby acknowledged and recognized as such. ©Copyright by Lifetrons Switzerland AG and Lifetrons International Limited. All rights are reserved.

"LIFETRONS" "LIETRONS" "LI" "LI-T" "LI-T" are trademarks or registered trademarks owned by Lifetrons Switzerland AG, Lifetrons International Ltd. or Li-Tek. All other brands and product names listed herein are trademarks or registered trademarks owned by their respective holders and are hereby acknowledged and recognized as such. ©Copyright by Lifetrons Switzerland AG, Lifetrons International Ltd. and Li-Tek. All rights reserved. Produced by: Li-Tek Electronics Technology Co., Ltd

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