

**Dot matrix RF Beauty detector**



All pictures shown are for illustration purpose only, actual products may vary.  
圖片僅供參考, 實物為主。  
圖片僅供參考, 實物為主。

RF-800A



**EN User Guide**

Thank you for purchasing this Lifetrans® product. This user guide outlines the operating procedures and safety precautions. Please read it carefully to make sure you are using your product properly and retain it for future reference.

**How to Use**

**A. RF & red light care**

1. Starting up: press the power button (3), buzzer "drops" as a reminder, dot matrix RF massage care open;
2. Please clean your skin before use. Apply skin cream (not included in the package) to the area you will be caring for. In order to achieve the best operation and obtain the best effect, we can also choose the special RF gel or cream (not included in the package). Press the power switch (3), start the equipment. The default gear setting is L (low) gear, press M function key (2) to adjust gear.
3. Hold the instrument to make sure the hands touch the conductive bar (5) on both sides and turn on to RF mode. The RF message head (4) must be fully contact with skin to turn on the red light function. Slide the device across the face with a little strength gently and slowly, massage over all areas of the face. (Suggestion: use for no more than 15 minutes each time and no more than 3 times a week.)
4. After 15 minutes of continuous use, the instrument will shut down automatically. If you want to shorten the using time, just hold the (3) 3 seconds, then shut down the instrument. Please clean the massage head with cotton pads after use.

**B. LED hint**

1. Low electricity prompt: LED lamp (1) flashing blue light.
2. When the power completely runs out, the product automatically shuts down and the buzzer "drops" three times.

When charging: LED light (1) always on.  
After the full of: LED lights (1) out.

**C. How to charge**

1. In the shutdown state, will (6) (Micro - USB plug) inserted into the port (6) (USB port) on the instrument;
2. Connect (8) (USB plug), to the USB ports which have power supply (for example: the AC converter with USB power supply), begin to charge;

**簡中 使用手冊**

感謝您購買 Lifetrans® 產品。本指南概述了操作方法和安全守則。請仔細閱讀以正確使用產品。請妥善保留本手冊以備將來參考之用。

**如何使用**

**A. 射頻&紅光護理**

1. 開機: 按一下電源鍵 (3), 蜂鳴器“滴”一聲提示, 點陣射頻按摩護理開啟;
2. 使用前請先清潔使用部位, 在即將護理的區域塗抹護膚霜 (包裝中不含)。為實現最佳操作, 獲得最佳效果, 亦可選用射頻專用凝膠或乳霜 (包裝中不含)。按一下電源開關 (3), 啟動儀器。默認檔位設置為L (低) 檔, 調節檔位則按M功能鍵 (2)。
3. 握住儀器, 確保手部接觸到兩側的導電棒 (5), 開啟進入射頻模式。射頻按摩頭 (4) 必須完全接觸皮膚, 以開啟紅光功能。稍稍用力, 輕柔緩慢地帶動儀器在臉部滑動, 按摩到所有部位。(建議: 每次使用時間不超過15分鐘, 每週使用不超過3次)
4. 持續使用15分鐘後, 儀器會自動關閉。如果想縮短使用時間, 只需按住 (3) 3秒, 即可關閉儀器。使用後請用化妝棉清潔按摩頭。

**B. LED提示**

1. 低電提示: LED (1) 發藍光閃爍。
2. 电量彻底用光時, 產品自動关机, 同時蜂鳴器“滴滴滴”三聲提示。

充電時: LED燈 (1) 常亮。  
充滿後: LED燈 (1) 熄滅。

**C. 如何充電**

1. 关机状态下, 將 (6) (Micro-USB插頭) 插入儀器上的端口 (6) (USB插口) ；
2. 連接 (8) (USB插頭), 到有電源供應的USB埠 (例如: 有USB電源供應的AC轉換器), 開始充電;

**繁中 使用手冊**

感謝您購買 Lifetrans® 產品。本指南概述了操作方法和安全守則。請仔細閱讀以正確使用產品。請妥善保留本手冊以備將來參考之用。

**如何使用**

**A. 射頻&紅光護理**

1. 開機: 按一下電源鍵 (3), 蜂鳴器“滴”一聲提示, 點陣射頻按摩護理開啟;
2. 使用前請先清潔使用部位, 在即將護理的區域塗抹護膚霜 (包裝中不含)。為實現最佳操作, 獲得最佳效果, 亦可選用射頻專用凝膠或乳霜 (包裝中不含)。按一下電源開關 (3), 啟動儀器。默認檔位設置為L (低) 檔, 調節檔位則按M功能鍵 (2)。
3. 握住儀器, 確保手部接觸到兩側的導電棒 (5), 開啟進入射頻模式。射頻按摩頭 (4) 必須完全接觸皮膚, 以開啟紅光功能。稍稍用力, 輕柔緩慢地帶動儀器在臉部滑動, 按摩到所有部位。(建議: 每次使用時間不超過15分鐘, 每週使用不超過3次)
4. 持續使用15分鐘後, 儀器會自動關閉。如果想縮短使用時間, 只需按住 (3) 3秒, 即可關閉儀器。使用後請用化妝棉清潔按摩頭。

**B. LED提示**

## Dot matrix RF Beauty detector



All pictures shown are for illustration purpose only, actual products may vary.  
图片仅供参考, 实物为准。  
RF-800A

## EN Safety Precautions

### Safety Precautions

Please follow the rules outlined in this guide closely to ensure your safety and flawless operation:

- Discontinue using the product immediately if any abnormal situation, allergic reaction or discomfort occurs during use.
- This product is a skin care device and is not made for medical purposes. If you are under medical treatment, please consult with your doctor before use.
- This product is not suitable to use on body parts that are embedded with silicone, metal or other foreign objects. Do not use directly on the eyes.
- Thoroughly clean the product after use.
- Do not use any kind of liquid to clean the product. Keep in mind that the surface of the product should only be wiped with a clean dry cloth.
- This product has a special lacquer - when cleaning, do not use any corrosive solution which may cause damage to the lacquer or appearance.
- Repair, removal or modification of this product is strictly prohibited. Improper handling may result in short circuits, product damage or even fire. Maintenance and warranty services will be void if unauthorized changes have been made to the product. Please contact an authorized dealer to obtain after-sales service if the product has a problem.
- Be cautious of any edges, uneven surfaces, metal parts, accessories or packaging to prevent possible injury or damage. Take note that this product should not be placed in your mouth.
- Do not use this product in locations with high temperatures or high humidity levels (e.g., bathrooms); ideally, the product should be used in a location free of excessive dust and not exposed to direct sunlight.
- The product is not waterproof. Do not insert foreign objects or liquids into the product; this could damage the product.
- Refrain from dropping the product or bumping it against other objects. Avoid damage, collision and excessive pressure by always handling the product with care.
- Do not place this product near strong magnetic or electric fields - this may cause product failure to occur.
- Do not use this product in hazardous locations.
- The product may explode if disposed in a fire.
- Due to hygienic reasons, do not share this product with other people.

### Children

Keep out of reach of children. Do not allow children to play the product or its accessories. They could hurt themselves or others and could accidentally damage the product. The product and its accessories may contain small parts it may detach and cause a choking hazard.

### Disposal of Electrical and Electronic Equipment

This symbol indicates that all electrical and electronic equipment included with this product should not be treated as household waste. Instead, it should be left at the appropriate collection point for the recycling of electronic equipment. By ensuring this product is disposed of correctly, you are helping to prevent negative consequences for the environment and human health. To recycle this product, please contact your local city office, environmental protection agency, household waste disposal service or where you purchased the product to learn more about appropriate recycling methods.

### Battery

It is recommended that you completely use the initial charge up, then recharge the product at least two more times to ensure the best possible performance.

A rechargeable battery can be charged and discharged over one hundred times. However, it will eventually wear out—this is not a defect and corresponds to normal everyday use. When the operating or standby time is noticeably shorter, it may be a sign that the battery has reached the end of its life span. Please do not make repairs, modifications, or installations without contacting an authorized Lifetrons representative.

### Limited Warranty

If within 12 months of purchase this product or any part thereof is proven to be defective by reason of faulty workmanship or defect in materials, we will repair or replace the same product free of charge of labor or materials on condition that:

1. The original invoice or sales receipt with the date of purchase is presented with the defective unit.
2. Lifetrons does not take any responsibility for lost information of any kind and will not reimburse you for any such loss.
3. The limited warranty is only valid if the original proof of purchase for this product, issued by a Lifetrons authorized dealer, specifying the date of purchase, is presented with the product to be repaired or replaced. Lifetrons reserves the right to refuse warranty service if this information has been removed or changed.
4. The warranty does not cover product failures which have been caused by the use of accessories or other peripheral devices which are not Lifetrons branded original accessories intended for use with the product. The warranty does not cover any damage caused by misuse, abuse, neglect or accidental breakage.
5. The warranty does not cover product failures caused by modifications, repairs or disassembly by any person not authorized by the manufacturer.
6. Any defective unit or part shall become property of the manufacturer.

### Additional Information

Register your Lifetrons Beauté product online at

[www.lifetrons.com](http://www.lifetrons.com).

Please visit our website [www.lifetrons.com](http://www.lifetrons.com) or contact our customer service center at [cs@lifetrons.com](mailto:cs@lifetrons.com) to learn more about our new products and special offers.

This user guide was originally in English and was subsequently translated into other languages. The accuracy of subsequent translations cannot be guaranteed. In case of conflict between the English version and another language version, the English version takes precedence.

Product specifications are subject to change at the own decision and discretion of Lifetrons Switzerland AG, Lifetrons International Ltd, or Li-Tek without prior notice. Any devices or accessories not listed under "What's Inside" are in reference to compatibility only and are not included. Product images are for reference purposes only, color may vary from the actual product.

"Lifetrons Beauté", "Lifetrons Beauté Switzerland", "Lifetrons Beauté Professional Skincare", "Lifetrons Switzerland", "Beauté Switzerland", "MyB" are trademarks or registered trademarks owned by Lifetrons Switzerland AG, Lifetrons International Ltd, or Li-Tek. All other brands and product names listed herein are trademarks or registered trademarks owned by their respective holders and are hereby acknowledged and recognized as such. ©Copyright 2018 by Lifetrons Switzerland AG, Lifetrons International Ltd, and Li-Tek. All rights reserved.

## 简中 安全守则

### 安全守则

为保安全和操作正确, 请先仔细阅读及遵守以下的安全守则:

- 如出现任何异常反应、过敏或不适, 请立即停止使用。
- 本产品只作护肤护理之用, 并非供医疗用途使用。如正接受任何治疗, 请先咨询医生。
- 本产品不适合于植入有硅胶、金属或其他外来物质的身体部分使用, 亦不可直接用于眼球上。
- 操作完毕后, 请保持产品的清洁。
- 请勿使用任何清洁剂或液体来清洁本产品的表面, 本产品只可使用干布清洁。
- 本产品有独特烤漆保护, 清洁时请勿使用任何腐蚀性溶液, 否则可能导致外表造成损坏。
- 本机严禁自行修理、拆卸或改装。不适当的处理可能导致电路短路、产品损坏甚至火灾。曾被擅自改动的产品, 均不保修。若产品产生问题, 请到指定售后单位获得售后服务。
- 任何边缘切割部位、不平滑的表面、金属零件、配件或包装均有可能因意外而引起流血或受伤, 请务必小心使用。不可将本产品置于口中。
- 请勿在高温或湿度高的地方(如浴室)使用本产品, 请在较少尘土的位置使用本产品, 不得暴露在诸如日照、火烤或类似的过热环境中。